

## Examination Rooms Item



*Excellence in Medical Technology*

### Introduction

WTK Manufacturing Sdn Bhd. (WTKM) is established in year 2005 under the roof of WTK Technologies (M) Sdn. Bhd. WTK Manufacturing Sdn Bhd. (WTKM) is dedicated to establishing a professional relationship with its clients and setting an example in hospital furniture system by providing today's hospital furniture system. WTKM will be the manufacturing arm to ensure the excellent quality of engineering work. We ensure you with a good market reputation and having good marketable records. Our Objective is to have a total customer satisfaction for our valued customers by provide the best solution of hospital furniture system by design & built, installation and excellent after sales service.

### Description

Technology and facilities - Our manufacturing processes are attributed to synergy of skilled personnel and advance technology. We are the first in our local industry to utilize state of art robotic welding systems, to ensure the highest quality of products. Today, almost all of our manufacturing processes are carried out with the help of computer aided machineries; starting from the design level with Auto-CAD to the fabrication level with computer programmed punching, CNC bending and auto-cutting machineries. In additional our business alliance and technological linkage with world class manufactures of enabling technologies have further enhanced our ability to deliver the superior quality products that complied with the international standard and safety requirements.

## WTK Manufacturing Sdn. Bhd.

6, Jalan BA/13, Kawasan Perusahaan Bukit Angkat

43000 Kajang, Selangor, Malaysia

Tel. : +603-87332880

Fax. : +603-87332535

E-mail : info@wtktechnologies.com.my



<http://www.wtktechnologies.com.my>

### PRODUCT SPECIFICATION AND PERFORMANCE

|   |              |   |
|---|--------------|---|
| <p><b>Examination Couch - Electrical System Back-up Battery</b></p> <ul style="list-style-type: none"><li>➤ Adjustable Backrest (40°) and Hi-Lo ( function through electric actuators</li><li>➤ Foot control pedal / foot switch.</li><li>➤ 100mm swivel castors with diagonal brake for mobility</li><li>➤ 50mm thick high density foam cushion padded with PVC upholstery padded.</li><li>➤ Optional Backrest with breathing hole with declinng(-20° )</li><li>➤ Approx. Dim(mm): 1830(L) x 600(W) x 475 - 995(H)</li></ul> | <b>E3001</b> |     |
| <p><b>Examination Table</b></p> <ul style="list-style-type: none"><li>➤ 2 section examination complete with foldable step stool</li><li>➤ Available EC1001(Without foot step)</li><li>➤ Approx. Dim(mm): 1830(L) x 600(W) x 760(H)</li></ul>  | <b>E1002</b> |     |
| <p><b>Examination Table</b></p> <ul style="list-style-type: none"><li>➤ 2 section examination</li><li>➤ C/w drawer and cabinet</li><li>➤ Optional with paper roll holder</li><li>➤ Approx. Dim(mm): 1830(L) x 600(W) x 760(H)</li></ul>   | <b>E2000</b> |   |
| <p><b>Examination table -Luxury Type</b></p> <ul style="list-style-type: none"><li>➤ Backrest adjustable</li><li>➤ C/w 3 drawers underneath</li><li>➤ C/w a mobile cabinet with 1 drawers</li><li>➤ C/w paper roll holder</li><li>➤ Approx. Dim(mm): 1830(L) x 600(W) x 760(H)</li></ul>  | <b>E2003</b> |  |
| <p><b>Examination Table –Gynae</b></p> <ul style="list-style-type: none"><li>➤ Backrest adjustable</li><li>➤ C/w stainless steel bowl</li><li>➤ C/w lithotomy stirrup</li><li>➤ Approx. Dim(mm): 1200(L) x 600(W) x 760(H)</li></ul>  | <b>E5000</b> |  |

TÜV Rheinland certified company for  
ISO 13485 & CE0197 Marking



**WTK Manufacturing Sdn. Bhd.**

6, Jalan BA/13, Kawasan Perusahaan Bukit Angkat

43000 Kajang, Selangor, Malaysia

Tel. : +603-87332880

Fax. : +603-87332535

E-mail : info@wtktechnologies.com.my



<http://www.wtktechnologies.com.my>

|  |   |   |
|--|---|---|
| <p><b>Gynae –Couch</b></p> <ul style="list-style-type: none"> <li>➤ Frame works steel with epoxy powder coated</li> <li>➤ Backrest adjustable by ratchet support</li> <li>➤ C/w knee cruth support adjustable rest</li> <li>➤ Foot rest fordable</li> <li>➤ C/w foot step</li> <li>➤ Approx. Dim(mm): 1200(L) x 600(W) x 760(H)</li> </ul> | <p><b>E5002</b></p>   |  |
| <p><b>Examination Light-LED</b></p> <ul style="list-style-type: none"> <li>➤ c/w a multi-flexible lamp arm.</li> <li>➤ LED bulb with 400 lumen; 25000 hours</li> <li>➤ 2700 kelvin</li> <li>➤ Voltage: 240V 50Hz</li> </ul>  | <p><b>M701<br/>LED</b></p>  |  |
| <p><b>X'Ray Viewer Box</b></p> <ul style="list-style-type: none"> <li>➤ Steel sheet constructed.</li> <li>➤ Wall-mounted/tabletop.</li> <li>➤ Finishing: Epoxy powder coated.</li> <li>➤ Voltage: 240V 50Hz. 32W lamp,</li> <li>➤ ON/OFF switch, 13A 3-pin plug,</li> <li>➤ shadow less perspex fascia</li> </ul>                          | <p><b>XRV101,<br/>XRV102,<br/>XRV103<br/>&amp;<br/>XRV104</b></p> |  |

TÜV Rheinland certified company for ISO 13485 & CE0197 Marking






# WTK Manufacturing Sdn. Bhd.

6, Jalan BA/13, Kawasan Perusahaan Bukit Angkat  
 43000 Kajang, Selangor, Malaysia  
 Tel. : +603-87332880  
 Fax. : +603-87332535  
 E-mail : info@wtktechnologies.com.my



<http://www.wtktechnologies.com.my>

|   |   |   |
|---|---|---|
| <p><b>LED X'Ray Illuminator</b></p> <ul style="list-style-type: none"> <li>➤ Thin X'Ray illumination ~ 25mm thickness</li> <li>➤ Separated control for different single viewer when it's more than single viewer</li> <li>➤ Super luminance is 5500cd/m2</li> <li>➤ Color temperature is over 9600k</li> <li>➤ Brightness can be controlled by knob or pressing key</li> <li>➤ The life of LED is 100000 hours</li> </ul> | <p><b>XRVL101,<br/>XRVL102,<br/>XRVL103<br/>&amp;<br/>XRVL104</b></p> |     |
| <p><b>Step Stool (Double Step)</b></p> <p><b>Material :</b> Mild Steel Framework<br/> <b>Finishing :</b> Epoxy powder coating<br/> <b>Dimension :</b> L500 x W230 x H330 mm<br/> <b>Step :</b> 1st H-180mm, 2nd H-330mm</p>   | <p><b>MST02</b></p>   |    |
| <p><b>Step Stool (Single Step)</b></p> <p><b>Material :</b> Mild Steel Framework<br/> <b>Finishing :</b> Epoxy powder coating<br/> <b>Dimension :</b> L500 x W230 x H180 mm<br/> <b>Step :</b> 1st H-180mm</p>  | <p><b>MST01</b></p>   |  |

TÜV Rheinland certified company for  
 ISO 13485 & CE0197 Marking

