

## Hospital Trolleys



### Introduction

WTK Manufacturing Sdn Bhd. (WTKM) is established in year 2005 under the roof of WTK Technologies (M) Sdn. Bhd. WTK Manufacturing Sdn Bhd. (WTKM) is dedicated to establishing a professional relationship with its clients and setting an example in hospital furniture system by providing today's hospital furniture system. WTKM will be the manufacturing arm to ensure the excellent quality of engineering work. We ensure you with a good market reputation and having good marketable records. Our Objective is to have a total customer satisfaction for our valued customers by provide the best solution of hospital furniture system by design & built, installation and excellent after sales service.

### Description

Technology and facilities - Our manufacturing processes are attributed to synergy of skilled personnel and advance technology. We are the first in our local industry to utilize state of art robotic welding systems, to ensure the highest quality of products. Today, almost all of our manufacturing processes are carried out with the help of computer aided machineries; starting from the design level with Auto-CAD to the fabrication level with computer programmed punching, CNC bending and auto-cutting machineries. In additional our business alliance and technological linkage with world class manufactures of enabling technologies have further enhanced our ability to deliver the superior quality products that complied with the international standard and safety requirements.

*Excellence in Medical Technology*

**PRODUCT SPECIFICATION AND PERFORMANCE**

<p><b>Medical Lab Dressing Table</b></p> <ul style="list-style-type: none"> <li>➤ Stainless Steel constructed.</li> <li>➤ Table top without guard rail</li> <li>➤ 2 shelves complete with 1 Drawer</li> <li>➤ Mobile on 75mm swivel castors</li> <li>➤ Dim(mm): 450W x 450D x 860H</li> <li>➤ <b>Material</b> : Stainless Steel 304 grade with hair line finishing.</li> </ul>	<p><b>T3201</b></p>	
<p><b>Medical Lab Dressing Table</b></p> <ul style="list-style-type: none"> <li>➤ Stainless Steel constructed</li> <li>➤ Table top without guard rail</li> <li>➤ 2 shelves complete with 2 Drawers</li> <li>➤ Mobile on 75mm swivel castors</li> <li>➤ Dim(mm):750W x 450D x 860H</li> <li>➤ <b>Material</b> : Stainless Steel 304 grade with hair line finishing.</li> </ul>	<p><b>T3202</b></p>	
<p><b>Medical Lab Dressing Table</b></p> <ul style="list-style-type: none"> <li>➤ Steel constructed with epoxy powder coated</li> <li>➤ Table top covered with stainless steel sheet and with 3 ways guard rail.</li> <li>➤ C/w 1 Drawer</li> <li>➤ mobile on 3" castors</li> <li>➤ Overall Dim(mm): 450(W) x 450(D) x 860(H) approx</li> </ul>	<p><b>T3001-E</b></p>	
<p><b>Medical Lab Dressing Table</b></p> <ul style="list-style-type: none"> <li>➤ Steel constructed with epoxy powder coated</li> <li>➤ Table top covered with stainless steel sheet and with 3 ways guard rail</li> <li>➤ C/w 2 Drawer</li> <li>➤ mobile on 3" castors</li> <li>➤ Overall Dim(mm): 750(W) x 450(D) x 860(H) approx</li> </ul>	<p><b>T3002-E</b></p>	

**WTK Manufacturing Sdn. Bhd.**

6, Jalan BA/13, Kawasan Perusahaan Bukit Angkat

43000 Kajang, Selangor, Malaysia

Tel. : +603-87332880

Fax. : +603-87332535

E-mail : info@wtktechnologies.com.my

<http://www.wtktechnologies.com.my>

<b>Medical Instrument Table</b> ➤ 2 shelves without guard rail Mobile on 75mm castors ➤ Overall Dim (mm): 450mm(L) x 450mm(W) x 860mm(H) approx.	<b>T2001</b>	
<b>Medical Instrument Table</b> ➤ 2 shelves without guard rail ➤ Mobile on 75mm castors ➤ Overall Dim (mm): 600mm (L) x 450mm(W) x 860mm(H) approx.	<b>T2002</b>	
<b>Medical Instrument Table</b> ➤ 2 shelves without guard rail ➤ Mobile on 75mm castors ➤ Overall Dim (mm): 750mm (L) x 450mm(W) x 860mmH approx.	<b>T2003</b>	
<b>Medical Instrument Table</b> ➤ 2 shelves without guard rail ➤ Mobile on 75mm castors ➤ Overall Dim (mm): 900mm (L) x 450mm(W) x 860mmH approx.	<b>T2004</b>	
<p>*** Optional with 3 ways guard rail with additional charges for each shelves</p>		
<b>Medical Instrument Table</b> ➤ Stainless steel constructed 304 grade ➤ 3 shelves with 3 ways bend up ➤ Mobile on 100mm castor ➤ Overall Dim(mm): 900mmL x 450mmW x 900mmH approx	<b>T2304</b>	

TÜV Rheinland certified company for ISO 13485 &amp; CE0197 Marking



<p><b>Mayo Stand Stainless Steel Constructed</b></p> <ul style="list-style-type: none"> <li>➤ chrome support rim for an s/steel tray</li> <li>➤ With adjustable height jam by screw knob.</li> <li>➤ Mobile on 50mm castors.</li> <li>➤ Tray size : - 490mmL x 320mmW</li> <li>➤ Overall Height: - 810 – 1230mmH approx</li> </ul>	<p><b>T4001</b></p>	
<p><b>COW-Lap Top</b></p> <ul style="list-style-type: none"> <li>➤ Steel construction with epoxy powder coated</li> <li>➤ Stand to accommodate Lap Top as the round trolley</li> <li>➤ Table top Hi-Lo adjustable for the comfortable of user</li> <li>➤ lap top compartment with key lock for the security of lap top</li> <li>➤ C/w sliding key board tray with built in mouse tray at side</li> <li>➤ C/w 4 in 1 extension cable and cable hook</li> <li>➤ Mobile on 4" twin wheels with 2 brakes</li> </ul>	<p><b>LTS01</b></p>	 
<p><b>COW-Desk Top</b></p> <ul style="list-style-type: none"> <li>➤ Steel construction with epoxy powder coated</li> <li>➤ Stand to accommodate desk top and CPU for paperless patient records in ICU / HDU unit</li> <li>➤ Table top Hi-Lo adjustable for the comfortable of user</li> <li>➤ Desk top table top with secure bar for the Monitor and base with CPU compartment to secure the CPU</li> <li>➤ C/w sliding key board tray with built in mouse tray at side</li> <li>➤ C/w single drawer below the sliding key board</li> <li>➤ C/w X-Ray films compartment \</li> <li>➤ C/w 4 in 1 extension cable and cable hook</li> <li>➤ Mobile on 4" twin wheels with 2 brakes</li> </ul>	<p><b>DTS01</b></p>	 

<p><b>Phebletomy Cart</b></p> <ul style="list-style-type: none"> <li>➤ Steel constructed with epoxy powder coated</li> <li>➤ Table top c/w lap top compartment with key lock</li> <li>➤ C/w sliding key board tray with built in mouse pad holder</li> <li>➤ C/w open compartment to accommodate printer</li> <li>➤ C/w 3 layer of drawer</li> <li>➤ side with sliding tray</li> <li>➤ C/w mask holder x 1 pcs</li> <li>➤ C/w Glove box holder x 1 pcs</li> <li>➤ C/w waste bin holder x 1 pcs</li> <li>➤ C/w sharp bin holder</li> <li>➤ C/w gel bottle holder x 1 pcs</li> <li>➤ Mobile on 4" twin wheels with 2 brakes</li> </ul>	<p><b>PBS01</b></p>	
<p><b>Drawer Cabinet For Oxygen Nebulizer</b></p> <ul style="list-style-type: none"> <li>➤ Steel constructed with epoxy powder coated</li> <li>➤ 2 layer of drawer different height</li> <li>➤ Bottom with cabinet with 2 doors</li> <li>➤ 4 leg with adjustable feet</li> </ul>	<p><b>CBO2N1</b></p>	
<p><b>ECG Cart</b></p> <ul style="list-style-type: none"> <li>➤ Steel constructed with epoxy powder coated</li> <li>➤ C/w with 1 drawer and and 1 shelve</li> <li>➤ Mobile on 4 x 50mm diameter castors</li> <li>➤ C/w with 1 hook</li> <li>➤ Approx Dim(mm): 400mm(L) x 350mm(W) x 850mm(H).</li> <li>➤ Optional Bottom with wiremesh basket</li> </ul>	<p><b>C6000</b></p>	
<p><b>Case Note Trolley</b></p> <ul style="list-style-type: none"> <li>➤ Mild steel plate construction</li> <li>➤ Scratch and chip resistant epoxy powder coating</li> <li>➤ Collapsible work table</li> <li>➤ Stainless steel rod and 15 PVC dividers for files</li> <li>➤ Mobile on 4 x 125mm dia. swivel castors with 2 brakes diagonally fitted</li> </ul>	<p><b>CCN01</b></p>	

## WTK Manufacturing Sdn. Bhd.

6, Jalan BA/13, Kawasan Perusahaan Bukit Angkat

43000 Kajang, Selangor, Malaysia

Tel. : +603-87332880

Fax. : +603-87332535

E-mail : info@wtktechnologies.com.my



<http://www.wtktechnologies.com.my>

### BHT Trolley

- Mild Steel constructed with epoxy powder coated
- Storage for folders/files at top 12 compartment
- Accommodate up to 60 files/folders (A4 size)
- 4 drawers
- X-Ray envelopes rack
- C/w pushing handle
- Mobility on 4 x 5" swivel castor with 2 brakes
- Collapsible working tray 22"L x 12"W
- C/w bumper
- Approx. Dim(mm) 720mm(L) x 560mm(W) x 950mm(H)

### CBHT1



TÜV Rheinland certified company for  
ISO 13485 & CE0197 Marking

