

## Hospital Carts



### Introduction

WTK Manufacturing Sdn Bhd. (WTKM) is established in year 2005 under the roof of WTK Technologies (M) Sdn. Bhd. WTK Manufacturing Sdn Bhd. (WTKM) is dedicated to establishing a professional relationship with its clients and setting an example in hospital furniture system by providing today's hospital furniture system. WTKM will be the manufacturing arm to ensure the excellent quality of engineering work. We ensure you with a good market reputation and having good marketable records. Our Objective is to have a total customer satisfaction for our valued customers by provide the best solution of hospital furniture system by design & built, installation and excellent after sales service.

### Description

Technology and facilities - Our manufacturing processes are attributed to synergy of skilled personnel and advance technology. We are the first in our local industry to utilize state of art robotic welding systems, to ensure the highest quality of products. Today, almost all of our manufacturing processes are carried out with the help of computer aided machineries; starting from the design level with Auto-CAD to the fabrication level with computer programmed punching, CNC bending and auto-cutting machineries. In additional our business alliance and technological linkage with world class manufactures of enabling technologies have further enhanced our ability to deliver the superior quality products that complied with the international standard and safety requirements.

*Excellence in Medical Technology*

**PRODUCT SPECIFICATION AND PERFORMANCE**

<p><b>Medication Cart</b></p> <ul style="list-style-type: none"> <li>➤ Construction-polymer mixed with aluminium and steel</li> <li>➤ Total of ABS bins 50 nos. (25 nos on front and 25 nos.at rear)</li> <li>➤ 5 layer ABS tray and each tray accommodate 5 ABS bins with divider inside</li> <li>➤ Bottom with big drawer with divider C/w 2 waste bins</li> <li>➤ C/w tillable hopper and top with bottle holder</li> <li>➤ ABS top with 3 ways guard rail</li> <li>➤ Mobile on 5" castors with 2 brakes</li> <li>➤ Overall Dim (mm): 850W x 520D x 950H approx</li> <li>➤ Optional with Overhead Bridge with 2 layer of dispenser and Rod with Hooks for model C3501</li> </ul>	<p>C3500</p>	
<p><b>Medication Cart 12 Large Bins</b></p> <ul style="list-style-type: none"> <li>➤ Stainless Steel construction for better infection control</li> <li>➤ Total 12 Bins large size (120W x 500D x 100H mm)</li> <li>➤ C/w 2 drawer at top and bottom</li> <li>➤ Mobile on 5" diameter castor with 2 brakes</li> <li>➤ C/w pushing handle</li> <li>➤ Overall dim(mm): 600W x 550D x 900H approx..</li> </ul>	<p>C3512</p>	
<p><b>Medical Emergency Cart -Polymer</b></p> <ul style="list-style-type: none"> <li>➤ Aluminium frame mild steel with epoxy and polymer constructed</li> <li>➤ 5 layer of drawers with divider inside</li> <li>➤ C/w defibrator board</li> <li>➤ C/w 2 waste bin</li> <li>➤ C/w tillable hopper and top with bottle holder</li> <li>➤ C/w IV pole with 4 hooks</li> <li>➤ C/w oxygen cylinder holder</li> <li>➤ C/ww CPR board</li> <li>➤ C/w sharp bin holder</li> <li>➤ Mobile on 5" castors with 2 brakes</li> <li>➤ Overall Dim(mm): 850W x 520D x 950H approx.</li> </ul>	<p>C3300</p>	

# WTK Manufacturing Sdn. Bhd.

6, Jalan BA/13, Kawasan Perusahaan Bukit Angkat  
43000 Kajang, Selangor, Malaysia  
Tel. : +603-87332880  
Fax. : +603-87332535  
E-mail : info@wtktechnologies.com.my



<http://www.wtktechnologies.com.my>

<p>Medical Emergency Cart</p> <ul style="list-style-type: none"> <li>➤ Steel constructed with epoxy powder coated</li> <li>➤ C/w 5 layer drawer</li> <li>➤ Table top with clear acrylic -see through window</li> <li>➤ 4 compartment bin on first 3 drawer</li> <li>➤ C/w defibrator tray</li> <li>➤ C/w oxygen cylinder holder</li> <li>➤ C/w utility tray beside oxygen cylinder holder</li> <li>➤ C/w catheter holder</li> <li>➤ C/w IV drip stand</li> <li>➤ C/w CPR Board</li> <li>➤ C/w all round bumper</li> <li>➤ Mobile on 5" castors with 2 brakes</li> <li>➤ Overall Dim(mm): 1000W x 460D x 970H approx.</li> </ul>	<p>C3100</p>	
<p>Medical Emergency Cart</p> <ul style="list-style-type: none"> <li>➤ Steel constructed with epoxy powder coated</li> <li>➤ C/w 4 layer drawer</li> <li>➤ Table top stainless steel constructed with or without epoxy powder coated</li> <li>➤ 4 compartment bin on first drawer</li> <li>➤ C/w collapsible defibrator tray</li> <li>➤ C/w oxygen cylinder holder</li> <li>➤ C/w IV drip stand</li> <li>➤ C/w CPR Board</li> <li>➤ C/w all round bumper</li> <li>➤ Mobile on 4" castors with 2 brakes</li> <li>➤ Overall Dim(mm): 1000W x 460D x 890H approx.</li> </ul>	<p>C3000</p>	
<p>Medical Injection Cart</p> <ul style="list-style-type: none"> <li>➤ Stainless steel constructed</li> <li>➤ Table top with 3 ways bend up</li> <li>➤ 4 layer of drawers equally divided</li> <li>➤ Each drawer with different size and number of compartments</li> <li>➤ C/w collapsible utility tray</li> <li>➤ C/w Pushing handle</li> <li>➤ C/w rotating corner bumper</li> <li>➤ Mobile on 5" castors with 2 brakes</li> <li>➤ Overall Dim(mm): 660W x 460D x 900H approx.</li> </ul>	<p>C3600</p>	

TÜV Rheinland certified company for ISO 13485 & CE0197 Marking



## WTK Manufacturing Sdn. Bhd.

6, Jalan BA/13, Kawasan Perusahaan Bukit Angkat

43000 Kajang, Selangor, Malaysia

Tel. : +603-87332880

Fax. : +603-87332535

E-mail : [info@wtktechnologies.com.my](mailto:info@wtktechnologies.com.my)



<http://www.wtktechnologies.com.my>

### File Trolley

- Steel constructed with epoxy powder coated
- Table top stainless steel constructed with 3 ways bend up
- C/w pushing handle
- C/w 10 arch file compartments
- C/w X-Ray films compartment
- C/w with utility compartment
- Mobile on 4 x 100mm diameter castors
- Overall Dim(mm): 900W x 460D x 900H approx.



TÜV Rheinland certified company for  
ISO 13485 & CE0197 Marking

