

Sleeper /Rehab / Rocking / Reclining Chair



Introduction

WTK Manufacturing Sdn Bhd. (WTKM) is established in year 2005 under the roof of WTK Technologies (M) Sdn. Bhd. WTK Manufacturing Sdn Bhd. (WTKM) is dedicated to establishing a professional relationship with its clients and setting an example in hospital furniture system by providing today's hospital furniture system. WTKM will be the manufacturing arm to ensure the excellent quality of engineering work. We ensure you with a good market reputation and having good marketable records. Our Objective is to have a total customer satisfaction for our valued customers by provide the best solution of hospital furniture system by design & built, installation and excellent after sales service.

Description

Technology and facilities - Our manufacturing processes are attributed to synergy of skilled personnel and advance technology. We are the first in our local industry to utilize state of art robotic welding systems, to ensure the highest quality of products. Today, almost all of our manufacturing processes are carried out with the help of computer aided machineries; starting from the design level with Auto-CAD to the fabrication level with computer programmed punching, CNC bending and auto-cutting machineries. In additional our business alliance and technological linkage with world class manufactures of enabling technologies have further enhanced our ability to deliver the superior quality products that complied with the international standard and safety requirements.

Excellence in Medical Technology

PRODUCT SPECIFICATION AND PERFORMANCE

| | | |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|--------------------------------------------------------------------------------------|
| <p>Medical Reclining / Rocking Chair</p> <ul style="list-style-type: none"> ➤ For use as Reclining / Relative Chair for patients on hemodialysis, chemotherapy, asthma, and recovery purpose ➤ With rocking function and base with 4 castors ➤ Back rest able to recline in to the most comfortable position and leg elevated when in reclining position by pulling the hand lever. ➤ Backrest 160° in position after recline. ➤ Material: ➤ Cushion with PVC padded ➤ Overall Dim(mm): ➤ 850 x 750D x 890H mm approx ➤ PVC colour selection: MOQ required | <p>CHR06</p> |  |
| <p>Reclining Chair</p> <ul style="list-style-type: none"> ➤ For use as Reclining / Relative Chair for patients on hemodialysis, chemotherapy, asthma, and recovery purpose. ➤ Reclining backrest at 160° in position after reclining and leg rest ➤ Material: Cushion with PVC padded c/w side pocket ➤ Overall Dim(mm): 800 W x 850 D x 1050 H approx. ➤ PVC colour selection: MOQ required | <p>CHR04</p> |  |
| <p>Geriatric Chair</p> <ul style="list-style-type: none"> ➤ Backrest adjustable ➤ PVC upholstery ➤ With detachable PU side table ➤ Mobile on 4 x 100mm diameter castors with 2 brakes ➤ PVC colour selection: MOQ required | <p>CHR07</p> |  |

| | | |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------|--------------------------------------------------------------------------------------|
| <p>Sleeper Chair</p> <ul style="list-style-type: none"> ➤ Sleeper Chair Cum Bed Converted from chair to bed ➤ Mobile on 75mm Twin wheel ➤ Frame work mild steel constructed with epoxy powder coated ➤ 3" cushion padded with PVC upholstery, with arm rest ➤ Seating Dim (mm): 720(W) x 900(D) x 1030H approx ➤ Sleeping Dim (mm): 1980(L) x 720(W) x 460(H) approx ➤ PVC colour selection: MOQ required | <p>CHSP02</p> |  |
| <p>Sleeper Chair / Accompany Chair</p> <ul style="list-style-type: none"> ➤ Sleeper Chair Cum Bed converted from chair to bed with Backrest Adjustable ➤ Backrest adjustable by ratchet support with various position ➤ Mobile on 75mm Twin wheel ➤ Frame work mild steel constructed with epoxy powder coated ➤ 3" cushion padded with PVC upholstery , with arm rest ➤ Seating Dim (mm): 720(W) x 900(D) x 1030H approx ➤ Sleeping Dim (mm): 1980(L) x 720(W) x 460(H) approx ➤ PVC colour selection: MOQ required | <p>CHSP04</p> |  |
| <p>Doctor Rotating Stool-Consultation room Without backrest</p> <ul style="list-style-type: none"> ➤ Seat area 350mm ➤ Hi-Lo adjustable 432-559 mm assisted by gas lift ➤ Aluminium base-chrome plated /aluminium base-epoxy mounted on 50mm castors ➤ Optional PU leg ➤ PVC colour selection: MOQ required. | <p>CHS01</p> |  |

| | | |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------|--------------------------------------------------------------------------------------|
| <p>ENT Doctor Chair</p> <ul style="list-style-type: none"> ➤ With semi backrest Hi-Lo adjustable 432-559 mm and assisted by gas lift ➤ Optional with footrest (CSH03) height adjustment 559-711mm ➤ PVC colour selection: MOQ required | <p>CHS02</p> |  |
| <p>Rotating High Lab Stool-Pneumatic</p> <ul style="list-style-type: none"> ➤ C/w Backrest Complete with footrest ➤ Mobile on PU / Aluminium leg mounted on 50mm castors. ➤ Hi-Lo adjustment by gas lift : 559-711mm ➤ Optional without footrest(CHS04) ➤ PVC Colour Selection: MOQ required | <p>CHS05</p> |  |
| <p>Surgeon Stool- foot Operated</p> <ul style="list-style-type: none"> ➤ C/w with backrest by single injection ➤ Aluminium base chassis ➤ Mobile on castors 50mm castors ➤ Seat Dim(mm) : 380mm ➤ Height: 510 - 660mm | <p>CHS06</p> |  |
| <p>Hospital Blood Taking Chair</p> <ul style="list-style-type: none"> ➤ Armrest adjustable and one side foldable ➤ Cushion padded with PVC upholstery ➤ Overall Dim(mm):650(650) x 520(D) x 900(H) | <p>CHBT01</p> |  |

WTK Manufacturing Sdn. Bhd.

6, Jalan BA/13, Kawasan Perusahaan Bukit Angkat

43000 Kajang, Selangor, Malaysia

Tel. : +603-87332880

Fax. : +603-87332535

E-mail : info@wtktechnologies.com.my



<http://www.wtktechnologies.com.my>

| | | |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <p>Hospital Blood Taking Chair</p> <ul style="list-style-type: none">➤ Backrest adjustable assisted by gas lift➤ Armrest Hi-Lo adjustable, tilt able and swing-able➤ Complete with hand grip for elderly or weak patient➤ Overall Dim(mm): 700(W) x 1000(D) x 1200(H) approx➤ Optional with leg rest and backrest manual adjustable by friction lock, mounted on castors➤ PVC colour selection: MOQ required | <p>CHBT02</p> |  Two photographs of the Hospital Blood Taking Chair (CHBT02). The left image shows the chair in its upright position, and the right image shows it reclined. The chair has a white metal frame, a blue padded seat and backrest, and adjustable armrests. |
| <p>Patient Chair / Visitor Chair</p> <ul style="list-style-type: none">➤ Wooden frame work with varnish coated➤ Good density cushion with PVC upholstery padded.➤ Overall Dim(mm): 655(W) x 610(D) x 1000(H) approx.➤ PVC colour selection: MOQ required | <p>CHPU01</p> |  A photograph of the Patient Chair / Visitor Chair (CHPU01). It is a simple, modern chair with a light-colored wooden frame and a white upholstered seat and backrest. |

TÜV Rheinland certified company for
ISO 13485 & CE0197 Marking

