

Patient Room Item



Introduction

WTK Manufacturing Sdn Bhd. (WTKM) is established in year 2005 under the roof of WTK Technologies (M) Sdn. Bhd. WTK Manufacturing Sdn Bhd. (WTKM) is dedicated to establishing a professional relationship with its clients and setting an example in hospital furniture system by providing today's hospital furniture system. WTKM will be the manufacturing arm to ensure the excellent quality of engineering work. We ensure you with a good market reputation and having good marketable records. Our Objective is to have a total customer satisfaction for our valued customers by provide the best solution of hospital furniture system by design & built, installation and excellent after sales service.

Description




Technology and facilities - Our manufacturing processes are attributed to synergy of skilled personnel and advance technology. We are the first in our local industry to utilize state of art robotic welding systems, to ensure the highest quality of products. Today, almost all of our manufacturing processes are carried out with the help of computer aided machineries; starting from the design level with Auto-CAD to the fabrication level with computer programmed punching, CNC bending and auto-cutting machineries. In additional our business alliance and technological linkage with world class manufactures of enabling technologies have further enhanced our ability to deliver the superior quality products that complied with the international standard and safety requirements.

Excellence in Medical Technology

PRODUCT SPECIFICATION AND PERFORMANCE

<p>Premier Over bed Table - Fingertip Height Adjustable</p> <ul style="list-style-type: none"> ➤ Table top laminated chipboard with ABS edging finishing . ➤ Table top Dim(mm): 900 x 400 approx ➤ Round edges cloumn laminated with stainless steel sheet with gold plated sheet with elegant look ➤ Hi-Lo adjustment with locking mecahnism to lock at variation of height just by 2 finger to operate it. ➤ Mobile on 50 mm Castors ➤ Dim(mm) : 900(L) x 400(W) x 760-1080(H) approx. 	<p>RX200</p>	
<p>Over bed Table - Fingertip Height Adjustable</p> <ul style="list-style-type: none"> ➤ Table top laminated chipboard with ABS edging finishing . ➤ Table top Dim(mm): 900 x 400 approx ➤ Hi-Lo adjustment with locking mecahnism to lock at variation of height just by 2 finger to operate it. ➤ Mobile on 50 mm Castors ➤ Overall Dim(mm) : 900(L) x 400(W) x 760-1080(H) approx. 	<p>RX202</p>	
<p>Over bed Table</p> <ul style="list-style-type: none"> ➤ Height Adjustable asissted by classic spring and screw knob ➤ Mobile on 50 mm Castors ➤ Table top laminated chipboard with ABS edging finishing ➤ Table top Dim(mm): 400 x 900 mm ➤ Overall Dim(mm) : 900(L) x 400(W) x 760-1080(H) approx. 	<p>RX201</p>	
<p>Overbed Table Fixed Height</p> <ul style="list-style-type: none"> ➤ Stainless steel constructed ➤ Table top with 3 ways bend up ➤ Mobile on 50mm castors ➤ Overall dim(mm): 1060(L) x 370(W) x 1000(H) approx ➤ Optional : with sharp bin holder 	<p>RX100-SS</p>	

<p>Hospital ABS Locker</p> <ul style="list-style-type: none"> ➤ C/w sliding tray x 1, drawer x 1, cabinet x 1, towel rail x 2 ➤ Overall dim(mm): 480W x 480D x 880H approx. ➤ Optional: Shoe Compartment and 50mm castor ➤ Colour selection: MOQ required 	<p>RX100</p>	
<p>Hospital Bedside Rack</p> <ul style="list-style-type: none"> ➤ Mild steel with epoxy powder coated ➤ Mobile on 50mm diameter castors ➤ C/w drawer x 1 and cabinet x 1 ➤ Overall Dim(mm): 460 x 400 x 812H approx. 	<p>RX101</p>	
<p>Hospital Bedside Rack</p> <ul style="list-style-type: none"> ➤ Mild steel with epoxy powder coated ➤ Table top with 3 ways guard rail ➤ C/w 1 drawer and 1 cabinet ➤ Overall Dim(mm): 403W x 520D x 835H approx. 	<p>RX102</p>	
<p>Hospital Bedside Rack</p> <ul style="list-style-type: none"> ➤ Mild steel with epoxy powder coated ➤ Mobile on 50mm diameter castors ➤ C/w drawer x 1 and cabinet x 1 ➤ Tap top dim(mm): 400(W) x 400(D) ➤ Overall Dim(mm): 400(W) x 560(D) x 812H approx. 	<p>RX103</p>	
<p>Bedside Locker</p> <ul style="list-style-type: none"> ➤ Chipboard laminated ➤ Complete with 1 drawer and 1 cabinet ➤ Rear with bottle compartment and towel rail ➤ Mobile on 4 x 50mm dia. Castor ➤ Overall Dim(mm): 400(W) x 445(D) x 810(H) approx 	<p>RX104</p>	

<p>Air Mattress (Transcell Type)</p> <ul style="list-style-type: none"> ➤ Heavy Duty ➤ C/w Pump <p>Pump Specification</p> <ul style="list-style-type: none"> ➤ Low energy consumption ➤ State-of-the-art pump with pressure control ➤ Light weight with hooks ➤ Easy to read, operate and clean <p>Mattress Specification</p> <ul style="list-style-type: none"> ➤ 4" replacement air cells ➤ Double layer designed pvc cell ➤ Low energy consumption ➤ Easy to clean cell surface ➤ 3 stages of alternating ➤ Water resistance and breathable mattress cover ➤ Material : pvc <p>Pump Technical Specification</p> <ul style="list-style-type: none"> ➤ Voltage : 110V/60Hz or 220V/50Hz ➤ Air out put : 4.5LPM ➤ Dim(mm) : 270(mm) x 140(W) x 80(H) <p>Mattress Technical Specification</p> <ul style="list-style-type: none"> ➤ Dim(mm) : 1800(L) x 840(W) x 110(H) 	<p>RXRM01</p>	
<p>Ripple Mattress (Bubbles Type)</p> <ul style="list-style-type: none"> ➤ Light Duty ➤ C/w pump ➤ Low Air Loss Mattress Specification ➤ Made of long lasting bubble textured PVC (EN 71, low toxic) ➤ Reusable, tear and accidental punctures resistance heavy gauge premium vinyl ➤ Height of cell : when inflated 2 1/2" <p>Pump Technical Specification</p> <ul style="list-style-type: none"> ➤ Voltage : 220V/50Hz ➤ Air out put : 5 – 6 liters ➤ Dim(mm) : 1982(L) x 834(W) (inflated) 	<p>RXRM02</p>	
<p>Hospital Mattress – Chemsafe / Mediflex</p> <ul style="list-style-type: none"> ➤ PU foam 32 Density ➤ Cover with water proof, anti-bacteria, anti-static and fire retardant ➤ Size: 1900mmL x 860mmW x 100mmH ➤ Colour: Beige / Blue 	<p>RXHM01</p>	

<p>Adult Pillow C/w Chemsafe Cover</p> <ul style="list-style-type: none"> ➤ Cover with Water-proof • Fire retardant • Anti-bacteria treated • Anti-static • Stain retardant ➤ Colour beige / Blue ➤ Overall dim(mm): 635(L) x 381(W) approx ➤ Also available for child RXHP02 <overall Dim(mm): 406(L) x 306(W) > approx 	<p>RXHP01</p>	
<p>Medical IV Drip Stand (2 hooks)</p> <ul style="list-style-type: none"> ➤ PU 5 Legged mounted on 50mm castors with SS Stand ➤ Stainless steel pole with 2 hook ➤ Pole adjustment: 1330 - 2300mm approx ➤ Also available pole with 4 hooks 	<p>MIVS01</p>	
<p>Medical IV Drip Stand (4 hooks)</p> <ul style="list-style-type: none"> ➤ Stand, pole and leg with stainless steel constructed. ➤ Pole with 4 hooks height adjustment 1330 - 2300mm approx.. ➤ Mobile on 3" castors 	<p>MIVS03</p>	
<p>Medical IV Drip Stand</p> <ul style="list-style-type: none"> ➤ PU 5 Legged mounted on 50mm castors with Epoxy coated stand ➤ Chrome plated pole with 2 hook ➤ Pole adjustment: 1330 - 2300mm approx 	<p>MIVSE1</p>	
<p>Hospital Ward Screen</p> <ul style="list-style-type: none"> ➤ 4 panels Screen with blue colour cloth ➤ Mobile on 50mm diameter twin castors ➤ Dim(mm): 2240L x 1850 ➤ Also available in 3 Panels or 2 Panels 	<p>MWS04</p>	