

Hospital Beds



Excellence in Medical Technology

Introduction

WTK Manufacturing Sdn Bhd. (WTKM) is established in year 2005 under the roof of WTK Technologies (M) Sdn. Bhd. WTK Manufacturing Sdn Bhd. (WTKM) is dedicated to establishing a professional relationship with its clients and setting an example in hospital furniture system by providing today's hospital furniture system. WTKM will be the manufacturing arm to ensure the excellent quality of engineering work. We ensure you with a good market reputation and having good marketable records. Our Objective is to have a total customer satisfaction for our valued customers by provide the best solution of hospital furniture system by design & built, installation and excellent after sales service.

Description

Technology and facilities - Our manufacturing processes are attributed to synergy of skilled personnel and advance technology. We are the first in our local industry to utilize state of art robotic welding systems, to ensure the highest quality of products. Today, almost all of our manufacturing processes are carried out with the help of computer aided machineries; starting from the design level with Auto-CAD to the fabrication level with computer programmed punching, CNC bending and auto-cutting machineries. In additional our business alliance and technological linkage with world class manufactures of enabling technologies have further enhanced our ability to deliver the superior quality products that complied with the international standard and safety requirements.

WTK Manufacturing Sdn. Bhd.

6, Jalan BA/13, Kawasan Perusahaan Bukit Angkat
43000 Kajang, Selangor, Malaysia
Tel. : +603-87332880
Fax. : +603-87332535
E-mail : info@wtktechnologies.com.my



<http://www.wtktechnologies.com.my>

PRODUCT SPECIFICATION AND PERFORMANCE

Hospitals Electric Bed-ICU / CCU (B6001)

Features :



- **HI-LO**, Backrest 75 deg., footrest 30 deg., electrically operated by embedded handset controller.
- Trendelenburg 15 deg. & reverse trend 15 deg., electrically operated by embedded handset controller.
- Quick release lever on bed frame, for emergency CPR
- Nurse Control unit at bed end for single or total lock of actuators.
- Night light on foot end of bed
- Backup Battery, for emergency power failure
- Split fold down safety side rail & Ventilated 4 panel mattress base.
- ABS easily removable head and foot board with lamination for easy mattress management.
- 125mm central locking twin wheel castors
- 2 IV pole mount attachment & roller bumpers at all corners to protect bed frame.
- SUS IV pole with double hook.
- Multi –purpose tray at base frame & urine bag holder at bed frame
- 100mm PU foam mattress with nylux cover.

Hospitals Electric Bed-Basic (B5001)

Features :

- **HI-LO**, Backrest 75 deg., footrest 45° elevation operated by electric Actuator.
- Back up battery
- Ventilated 4-panel mattress base
- ABS moulded head and foot board with Maica lamination.
- 2 diagonally fitted brakes on swivel castors
- 1 unit of IV pole. Model IVP1 come with IV pole mount at bed frame
- Roller bumpers at all corners, to protect body frame
- Collapsible Epoxy finish safety rails with spring loaded lock.
- 100mm PU foam Mattress with “MEDISOFT” Nylux covering.



Hospital Hydraulic Bed (B4001)

Features :



- Backrest (75°) elevation by Gas Lift.
- Footrest (45°) elevation by Manual Hand Crank
- Hi-Lo elevation operated by hydraulic foot pedal
- 2-way Tilting – Trendelenburg (12°) & Reverse Trend (6°) by gas activation.
- Ventilated 4-panel mattress base
- ABS-moulded head and foot board, with Maica lamination (Made in Malaysia)
- 125mm swivel castor with Central Locking.
- Roller bumpers at all corners, to protect body frame.
- Collapsible Chromed finish safety side rails with spring-loaded lock

TÜV Rheinland certified company for
ISO 13485 & CE0197 Marking



WTK Manufacturing Sdn. Bhd.

6, Jalan BA/13, Kawasan Perusahaan Bukit Angkat

43000 Kajang, Selangor, Malaysia

Tel. : +603-87332880

Fax. : +603-87332535

E-mail : info@wtktechnologies.com.my



<http://www.wtktechnologies.com.my>

Hospital Ultra Low Bed (B5002)

Features :

- HI-LO, backrest (75°) & footrest (45°) elevation operated by Electric Actuator.
- Back-up battery.
- Ventilated 4-panel mattress base
- ABS-moulded head and foot board, with Maica lamination
- 2 diagonally fitted brakes on swivel castors
- 1 unit of IV pole Model: IVP1 c/w IV pole mount at bed frame
- Roller bumpers at all corners, to protect body frame
- Collapsible Epoxy-finish safety side rails with spring-loaded lock
- 100mm PU foam Mattress with "MEDISOFT" Nylex covering.



Hospital Manual Crank Bed (B3001)

Features :

- Backrest, legrest and hi-Lo by manual cranking
- Comes with telescopic side rail, mobile on castors with 2 brakes and also comes with hospital grade mattress.



Maternity/ Delivery Beds Electric System (LDR30)

Features :

Backrest elevation angle up to 70° and Hi-Lo adjustment by electric actuators, operated by foot switch.



Trendelenburg position from 0 ~ 12°. 125mm central locking swivel castors. This unit comes with hand grip, stainless steel bowl, knee crutch and IV pole.



Convertible length of stainless steel leg frame from bed to delivery unit.



Adjustable backrest to prepare for delivery position.



Trendelenburg position from 0 – 12 degree.



Adjustable leg rest and foot stirrup and hand grip



Rotate down safety side rails



Stainless steel fluid collection bowl

TÜV Rheinland certified company for ISO 13485 & CE0197 Marking



TÜV Rheinland®
CERT
EN ISO 13485



WTK Manufacturing Sdn. Bhd.

6, Jalan BA/13, Kawasan Perusahaan Bukit Angkat

43000 Kajang, Selangor, Malaysia

Tel. : +603-87332880

Fax. : +603-87332535

E-mail : info@wtktechnologies.com.my



<http://www.wtktechnologies.com.my>

Maternity/ Delivery Beds Hydraulic-System (LDR20)

Features :

Backrest elevation angle up to 70° by gas lift and Hi-Lo adjustment by hydraulic operated by foot pedal on both sides .

Gas spring manual trendelenburg (tilting) : 0 ~ 12°. This unit is come with hand grip, stainless steel bowl, knee crutch and IV pole.



****Note : Both units come with 100mm thickness of PU mattress**

TÜV Rheinland certified company for
ISO 13485 & CE0197 Marking

